



659 Andover Street
Lawrence, MA 01843-1032, USA
Tel: 1.978.687.9200 Fax: 1.978.687.9216

Certificate of Quality

Product Name: Charm ROSA Tetracycline Test for Milk at U.S. Tolerance

Lot Number: TET-SL-112

Expiration Date: October 2023

Order Codes: LF-TET-SL

This product was quality tested and met Charm Sciences' quality control criteria for release for sale including performance and appearance specifications.

Name: Sam Atoui
Title: Quality Assurance Documentation Manager
Date: 03/15/2023

CHARM TET-SL – U.S. Shipping Record

CHARM TET-SL MUST BE RECEIVED WITHIN 72 HOURS IF SHIPPED NON REFRIGERATED

CUSTOMER USE ONLY:

SHIP DATE (from packing slip)	RECEIVED DATE (within 72 hours of SHIP DATE)	STORED REFRIGERATED	INITIALS
_____	_____	_____	_____

Nothing works like a Charm



659 Andover Street
Lawrence, MA 01843-1032, USA
Tel: 1.978.687.9200 Fax: 1.978.687.9216

Certificate of Quality

Product Name: Tetracycline Dilution Buffer

Lot Number: TDB-005

Expiration Date: December 2023

Order Code: BUF-TDB

This product was quality tested and met Charm Sciences' quality control criteria for release for sale including performance and appearance specifications.

A handwritten signature in black ink, appearing to read "Sam Atoui", is written above a horizontal line.

Name: Sam Atoui

Title: Quality Assurance Documentation Manager

Date: 08/31/2022

TYPICAL COUNTS for
Charm SL Tetracycline Test

See Operator's Manual for use and application.

Kit Code: LF-TET-SL

Lot Number: TET-SL-112

Charm ROSA Reader

Negative Control	-2854	3434	MRLBLTET/TET Positive Tablet (300 ppb Oxytetracycline)
	-2465	3530	
	-2478	3142	
	-3505	2891	
	-2938	2452	
	-2432	2348	
	-2980	2575	
	-3047	2641	
	-2593	1913	
	-3031	2709	

Note: Typical counts are provided as a performance guide only.
Operators' daily control numbers may be compared to the supplied
information as an operator, equipment and analyzer calibration check. Daily
control data should typically test in the range of data points provided.
The control point for the ROSA SL TetracyclineTest is always zero.