

**TYPICAL COUNTS for**  
**Charm QUAD2 Test for Macrolides**

See Operator's Manual for use and application.

Kit Code: LF-QUAD2

**Lot Number:      QUAD2-007AA**

Charm ROSA EZ Reader

Negative Control	-3587	3129	QUAD2 Positive Control
	-3591	2909	200 ppb Spiramycin,
	-2663	2656	80 ppb Erythromycin,
	-3222	2856	300 ppb Lincomycin, and
	-3042	3237	200 ppb Tylosin
	-3436	2903	
	-3092	2914	
	-3411	3178	
	-3732	3156	
	-3849	2842	

Note: Typical counts are provided as a performance guide only.

Operators' daily control numbers may be compared to the supplied information as an operator, equipment and analyzer calibration check. Daily control data should typically test in the range of data points provided.

The control point for the QUAD2 Test is always zero.

02/07/20