

### **Certificate of Quality**

Product Name: Charm ROSA Tetracycline Test for Milk at U.S. Tolerance

Lot Number: TET-SL-111

Expiration Date: September 2023

Order Codes: LF-TET-SL

This product was quality tested and met Charm Sciences' quality control criteria for release for sale including performance and appearance specifications.

Name: Sam Atoui

Title: Quality Assurance Documentation Manager

Date: 02/07/2023

# **CHARM TET-SL – U.S. Shipping Record**

\*CHARM TET-SL MUST BE RECEIVED WITHIN 72 HOURS IF SHIPPED NON REFRIGERATED\*

CUSTOMER USE ONLY:

SHIP DATE (from packing slip)

RECEIVED DATE (within 72 hours of SHIP DATE)

STORED REFRIGERATED

**INITIALS** 



## **Certificate of Quality**

Product Name: Tetracycline Dilution Buffer

Lot Number: TDB-005M

Expiration Date: December 2023

Order Code: BUF-TDB

This product was quality tested and met Charm Sciences' quality control criteria for release for sale including performance and appearance specifications.

Name: Sam Atoui

Title: Quality Assurance Documentation Manager

Date: 11/03/2022

#### **TYPICAL COUNTS for**

### **Charm SL Tetracycline Test**

See Operator's Manual for use and application.

Kit Code: LF-TET-SL

Lot Number: TET-SL-111

Charm ROSA Reader

Negative Control	-2854	3434	MRLBLTET/TET Positive Tablet
	-2465	3530	(300 ppb Oxytetracycline)
	-2478	3142	
	-3505	2891	
	-2938	2452	
	-2432	2348	
	-2980	2575	
	-3047	2641	
	-2593	1913	
	-3031	2709	

Note: Typical counts are provided as a performance guide only.

Operators' daily control numbers may be compared to the supplied information as an operator, equipment and analyzer calibration check. Daily control data should typically test in the range of data points provided.

The control point for the ROSA SL TetracyclineTest is always zero.